

Head-to-head comparison of the fully automated Elecsys pTau217 plasma assay and the Lumipulse pTau217 plasma assay

Robert Pernecky,^{1–5} Laura Stoeckl,⁶ Margherita Carboni,⁷ Sayuri Hortsch,⁶ Christina Rabe,⁸ Alexander Jethwa,⁶ Tobias Bittner^{8,9}

¹University Hospital of Munich, Munich, Germany; ²German Center for Neurodegenerative Diseases, Munich, Germany; ³Munich Cluster for Systems Neurology, Munich, Germany; ⁴Ageing Epidemiology Research Unit, London, UK; ⁵Sheffield Institute for Translational Neuroscience, Sheffield, UK; ⁶Roche Diagnostics GmbH, Penzberg, Germany; ⁷Roche Diagnostics International Ltd, Rotkreuz, Switzerland; ⁸Genentech Inc., South San Francisco, USA; ⁹F. Hoffmann-La Roche Ltd, Basel, Switzerland

Poster no. LP069

Introduction

- Blood tests are a non-invasive and cost-effective method for detecting Alzheimer's disease pathology, including amyloid deposits.¹
- Plasma phosphorylated tau 217 (pTau217) is a biomarker that has demonstrated high clinical utility in detecting amyloid pathology compared with other biomarkers.^{2,3}
 - The Elecsys® Phospho-Tau 217 plasma prototype immunoassay (Roche Diagnostics International Ltd, Rotkreuz, Switzerland) quantitatively determines plasma pTau217 levels to identify amyloid pathology and has received United States Food and Drug Administration Breakthrough Device Designation.⁴
 - The Lumipulse G pTau 217 plasma assay is a research use only assay from Fujirebio.⁵

Objective

- To perform a head-to-head method comparison study of the Elecsys and Lumipulse pTau217 assays.

Methods

- This head-to-head, blinded, method comparison study compared plasma pTau217 levels measured using the Elecsys and Lumipulse pTau217 assays in samples from a selected subset of patients from the amyloid screening population of the CREAD2 study (NCT03114657).⁶
- Amyloid status was determined by amyloid-positron emission tomography (PET) imaging with visual read or, if not available, by the cerebrospinal fluid (CSF) pTau181/β-amyloid(1–42) ratio, measured with the Elecsys Phospho-Tau (181P) CSF and Elecsys β-Amyloid(1–42) CSF immunoassays (Roche Diagnostics International Ltd, Rotkreuz, Switzerland).
- The receiver operating characteristic (ROC) area under the curve (AUC) with respect to amyloid status was calculated.
- A dual cutoff approach assessed the cutoffs at sensitivity and specificity of 90%. Positive percent agreement (PPA), negative percent agreement (NPA), positive predictive values (PPV) and negative predictive values (NPV), were calculated.
- The Pearson correlation coefficient between results from the two pTau217 assays was calculated.

Results

- This analysis included samples from 264 patients (80 [30.3%] were amyloid-negative and 184 [69.7%] were amyloid-positive; **Table 1**).

Table 1. Patient demographics and baseline characteristics by amyloid status.

	Amyloid-negative (n=80)	Amyloid-positive (n=184)	All (N=264)
Mean age (range), years	72.4 (50.0–85.0)	70.7 (53.0–85.0)	71.3 (50.0–85.0)
Female, n (%)	40 (50.0)	106 (57.6)	146 (55.3)
Race, n (%)			
White	69 (86.3)	166 (90.2)	235 (89.0)
Other*	11 (13.8)	18 (9.8)	29 (11.0)
Mean BMI (range), kg/m²	27.1 (20.0–36.8)	25.4 (15.8–42.4)	25.5 (15.8–42.4)
APOE E4 alleles, n (%)			
No copies	52 (65.0)	58 (31.5)	110 (41.7)
One copy	22 (27.5)	98 (53.3)	120 (45.5)
Two copies	6 (7.5)	28 (15.2)	34 (12.9)
Mean MMSE score (range)	25.1 (22.0–29.0)	24.5 (22.0–30.0)	24.7 (22.0–30.0)
PET status, n (%)			
Negative	61 (76.3)	0 (0.0)	61 (23.1)
Positive	0 (0.0)	147 (79.9)	147 (55.7)
Missing	19 (23.8)	37 (20.1)	56 (21.2)

*Other races included American Indian or Alaska Native, Asian, Black or African American, multiple races or unknown. APOE E4, Apolipoprotein E4; BMI, body mass index; MMSE, Mini-Mental State Examination.

- For the Elecsys pTau217 assay, six (2.3%) measurements were below the corresponding limits of quantification compared with one (0.4%) measurement for the Lumipulse pTau217 assay, demonstrating sensitivity to measure pTau217 within the cohort of samples.
- For the Elecsys pTau217 assay, there were no measurements that exceeded the assay measuring range compared with two (0.8%) measurements for the Lumipulse pTau217 assay.

Acknowledgements and Disclosures

The authors would like to acknowledge Judith Schmitz (Chrestos Concept GmbH & Co. KG, Essen, Germany) for contributing to the data analysis for this study. R.P. has received honoraria for advisory boards and speaker engagements from Roche, Eisai, Eli Lilly, Biogen, Janssen-Cilag, AstraZeneca, Schwabe, Grifols, Novo Nordisk, GSK, Bristol Myers Squibb and Tabak. L.S., S.H. and A.J. are full-time employees of Roche Diagnostics GmbH. M.C. is a full-time employee of Roche Diagnostics International Ltd and a shareholder in Roche. C.R. is a full-time employee of Genentech Inc., a member of the Roche Group and shareholder in F. Hoffmann-La Roche AG. T.B. is a full-time employee of F. Hoffmann-La Roche Ltd and Genentech Inc., a member of the Roche Group and shareholder in F. Hoffmann-La Roche AG. This study was funded by Roche Diagnostics International Ltd (Rotkreuz, Switzerland). Third-party medical writing assistance for the development of this poster, under the direction of the authors, was provided by Tiffany Blythe, BSc, of Ashfield MedComms, an Intizo company (London, UK), and funded by Roche Diagnostics International Ltd (Rotkreuz, Switzerland). ELECSYS is a trademark of Roche. All other product names and trademarks are the property of their respective owners. The Elecsys Phospho-Tau 217 prototype plasma immunoassay is not currently approved for clinical use. The Elecsys Phospho-Tau (181P) CSF and Elecsys β-Amyloid(1–42) CSF immunoassays are approved for clinical use.

Results (cont.)

- Plasma pTau217 levels between the Elecsys and Lumipulse pTau217 assays were positively correlated (Pearson's $r=0.883$; ten outliers identified by Generalized Extreme Studentized Deviate were removed; **Figure 1**).
- The AUC by amyloid status for pTau217 was 0.907 (95% confidence interval [CI]: 0.858–0.957) with the Elecsys pTau217 assay and 0.862 (95% CI: 0.805–0.920) with Lumipulse (**Figure 2**).

Figure 1. Correlation of plasma pTau217 levels measured using Elecsys and Lumipulse pTau217 assays by amyloid status.






Figure 2. ROC curves with AUC for pTau217 with respect to amyloid status* measured using the Elecsys and Lumipulse pTau217 assays.



*Sensitivity/specificity were calculated with respect to amyloid status, determined by amyloid-PET visual read, if available, or otherwise using the Elecsys CSF pTau181/β-amyloid(1–42) ratio.

- Both the Elecsys and Lumipulse pTau217 assays demonstrated good discrimination of amyloid-negative and amyloid-positive samples with the 90% sensitivity and 90% specificity cutoffs (**Figure 3**).
 - For the Elecsys pTau217 assay, 10.2% of samples were in the intermediate zone between the two cutoffs compared with 28.0% with the Lumipulse pTau217 assay.

Figure 3. Dual cutoff approach of 90% sensitivity (orange line) and 90% specificity (red line) for the (A) Elecsys pTau217 assay and (B) Lumipulse pTau217 assay*.



*Each box represents the interquartile range (25th–75th percentile) with the horizontal bold line showing the median value and whiskers representing the upper and lower boundaries for the extreme limits. The orange line represents the 90% sensitivity cutoff, and the red line represents the 90% specificity cutoff.

Conclusions

- The Elecsys pTau217 assay detects amyloid pathology with high accuracy and is highly correlated with the Lumipulse pTau217 assay.
- The fully automated Elecsys and Lumipulse pTau217 assays have clinical potential for detecting amyloid pathology. However, the assays need to be fully validated and approved by regulatory bodies before clinical use.

References

- Hampel H, et al. Nat Aging 2022;2:692–703.
- Palmqvist S, et al. JAMA 2020;324:772–81.
- Thijssen EH, et al. Lancet Neurol 2021;20:739–52.
- Roche. Roche granted FDA Breakthrough Device Designation for blood test to support earlier Alzheimer's disease diagnosis; 2024. Available at: <https://www.roche.com/media/releases/med-cor-2024-04-11>. Accessed 24 September 2024.
- Fujirebio. Lumipulse G pTau 217 Plasma. Available at: <https://www.fujirebio.com/en/products-solutions/lumipulse-g-ptau-217-plasma>. Accessed 24 September 2024.
- ClinicalTrials.gov. NCT03114657. Available at: <https://www.clinicaltrials.gov/study/NCT03114657>. Accessed 24 September 2024.

Access this poster using the quick response (QR) code or by visiting:

<https://ter.li/5wjudh>

